

LACTOSCAN TBS

Thermostatic device / Dry incubator

Operation manual

1. INTRODUCTION

LACTOSCAN TBS is very simple to use thermostatic device for incubation of BetaStar® Rapid Test and BetaStar® combo Rapid Test. The users only have to choose the desired incubation mode.

For safe and dependable operation of this thermostatic device, please comply with the following **safety precautions**:

- Verify that the input voltage printed on the AC Adapter and the plug type matches the local AC power supply.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Place the incubator on a firm and level surface.
- Disconnect the thermostatic device from the power supply when cleaning it.
- Follow the instructions in this manual when operating with the thermostatic device.
- Service should be performed by authorized personnel only.

2. Unpacking:

Unpack and verify that the following components have been included:

- LACTOSCAN TBS;
- Switching AC adapter 12V or 24V;
- Power cord;
- Car/truck lighter connecting cord;
- Additional cover without holes – to be used when the ambient temperature is below 10 °C;
- Instruction manual.

Save the packaging materials. This packaging ensures the best possible protection for the storage or transport of the product.

3. Mechanical parameters:

W x L x H: 90 x 157 x 69 mm

Weight: 0,780 kg

Stainless steel cover box




4. Electrical parameters:

LACTOSCAN TBS can be supplied with 12V or 24V power supply. On customers' request the device is shipped with either 12V or 24V switching adapter.

Switching adapter details:

- **12V switching adapter:**
Input: 100-240V ~ 50/60 Hz
Output: 12V $\ddot{=}$ 3A
- **24V switching adapter:**
Input: 100-240V ~ 50/60 Hz
Output: 24V $\ddot{=}$ 3.75A

5. Overview of display indications and keyboard functions:



- Timer 1 (green) – shows the status of Timer 1 (also in the programming mode).
- Timer 2 (green) – shows the status of Timer 2 (also in programming mode).
- Heater (red) – shows the status of the heater (on-off).
- Button SET – serves to adjust the temperature / time.
- Button START – starts Timer 1 and Timer 2.
- Three digits indication (display) - shows the temperature / time.
- Decimal point – blinks when measuring the temperature.

6. Operations:

6.1. Factory settings:

The thermostatic device is factory set for the following temperature and time ranges in Mode 1 and Mode 2:


	Mode 1 BetaStar® Rapid Test	Mode 2 BetaStar® combo Rapid Test
Temperature	47.5 °C	47.5 °C
Timer 1	3 minutes	2 minutes
Timer 2	2 minutes	3 minutes

7. Switching modes of LACTOSCAN TBS and starting the incubation:

7.1. Selecting Mode 1:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press and hold the button **SET**.
- Switch on the thermostatic device.
- After 3 seconds the device beeps and on the display is written the number "1", the software version "In6" and 475.
- Release the button SET.
- The device restarts and on the display is the current temperature.


See photos 1-4 below:



7.2. Selecting Mode 2:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press and hold the button **START**.
- Switch on the thermostatic device.
- After 3 seconds the device beeps and on the display is written the number "2", the software version "In6" and 475.
- Release the button START.
- The device restarts and on the display is the current temperature.

See photos 1-4 below:




After selecting the desired incubation mode, the thermostatic device is ready to be used.

8. Starting incubation:

- After selecting the desired incubation mode, wait till the thermostatic device reaches the set temperature for the mode.
- Once the temperature is reached, press **START** button to start running Timer 1.
- When Timer 1 finishes, the device beeps, on the display is written number "1" and the user has to press the **START** button to start Timer 2.
- When Timer 2 finishes, the device beeps, on the display is written number "2" and the user has to press the **START** button to stop the mode.

See photos 1-11 below:



9. Programming:


For programming or reprogramming, the user must follow the step below:

9.1. Setting parameters of Mode 1:

9.1.1. Choosing Mode 1:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press and hold the button **SET** and keep it pressed.
- Switch on the thermostatic device.
- After 3 seconds the device beeps and on the display is written the number "1", the software version "In6" and 475.
- Release the button SET.
- The device restarts and on the display is the current temperature.

See photos 1-4 below:




9.1.2. Setting temperature for Mode 1:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press the buttons **SET** and **START** simultaneously and keep them pressed.
- Switch on the thermostatic device.
- After 2 seconds, on the display is written inscription SET accompanied by a double beep.
- Release both buttons.
- The display shows 30.0 – this is the initial temperature.

- By using the **SET** button, the temperature can be changed upward to 80 degrees.
- Press the button **START** to save the adjusted temperature.

See photos 1-7 below:



9.1.3. Setting time for Timer 1 of Mode 1:

- After you have saved the adjusted temperature, on the display is written 001. Use the button **SET** to adjust the desired time for Timer 1. By pressing the button **SET** the time interval is changed to 1 minute.
- Press the button **START** to save the adjusted time.

See photos 1-3 below:



9.1.4. Setting time for Timer 2 of Mode 1:

- After you have saved the adjusted time for Timer 1, on the display is written 001. Use the button **SET** to adjust the desired time for Timer 2. By pressing the button **SET** the time interval is changed to 1 minute.
- Press the button **START** to save the adjusted time.

See photos 1-3 below:



The device beeps and shows the current temperature. Now the temperature and time settings for mode 1 have been saved.


9.2. Setting parameters of Mode 2:

9.2.1. Choosing Mode 2:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press and hold the button **START** and keep it pressed.
- Switch on the thermostatic device.

- After 3 seconds, the device beeps and on the display is written number “2”, the software version “In6” and 475.
- Release the button START.
- The device restarts and on the display is the current temperature.

See photos 1-4 below:



1.

2.


3.

4.

9.2.2. Setting temperature for Mode 2:

- Switch off the thermostatic device from the ON / OFF Power button.
- Press the buttons SET and START simultaneously and keep them pressed.
- Switch on the thermostatic device.
- After 2 seconds, on the display is written inscription SET accompanied by a double beep.
- Release both buttons.
- The display shows 30.0 – this is the initial temperature.
- By using the SET button, the temperature can be changed upward to 80 degrees.
- Press the button START to save the adjusted temperature.

See photos 1-7 below:



1.

2.

3.

4.



5.

6.

7.

9.2.3. Setting time for Timer 1 of Mode 2:

- After you have saved the adjusted temperature, on the display is written 001. Use the button SET to adjust the desired time for Timer 1. By pressing the button SET the time interval is changed to 1 minute.
- Press the button START to save the adjusted time.

See photos 1-3 below:



1.

2.

3.

9.2.4. Setting time for Timer 2 of Mode 2:

- After you have saved the adjusted time for Timer 1, on the display is written 001. Use the button SET to adjust the desired time for Timer 2. By pressing the button SET the time interval is changed to 1 minute.
- Press the button START to save the adjusted time.

See photos 1-3 below:



1.

2.

3.

The device beeps and shows the current temperature. Now the temperature and time settings for mode 2 have been saved.

10. Connecting LACTOSCAN TBS towards 12V car lighter or 24V truck lighter:

The user can connect the device towards a car lighter voltage 12V or towards a truck lighter voltage 24V by using the car/truck lighter connecting cord.

11. Cleaning:



ATTENTION!

In order not to damage any part of the thermostatic device, please follow the cleaning instructions below.

Clean the **metal cover box** of the thermostatic device with ethyl alcohol and soft cloth.

Clean the **plastic parts and covers** of the thermostatic device with dishwashing liquid, water and soft cloth.

LACTOSCAN Thermostatic device Data sheet

Model:	LACTOSCAN TBS
Number of wells:	13
Compatible test/tube:	Beta Star® Rapid Test and Beta Star® combo Rapid Test
Temperature range:	30°C to 80°C ($\pm 0.5^\circ\text{C}$)
Working modes:	1 st mode: <ul style="list-style-type: none"> • Temperature: 47.5 °C • 1st time interval: 3 minutes • 2nd time interval: 2 minutes
	2 nd mode: <ul style="list-style-type: none"> • Temperature: 47.5 °C • 1st time interval: 2 minutes • 2nd time interval: 3 minutes
Environmental condition:	Ambient temperature has to be 5°C below the desired set temperature
	Operating temperature: 0 to 40°C
Electrical parameters:	Switching adapter: Input 100 – 240 V ~ 1.6 Amax, 50 ~ 60 Hz
	Power consumption: 25W (12V), 100W (24V) until the necessary temperature is reached, after that 1/10 th
Dimensions (W x L x H):	90 x 157 x 69 mm
Weight:	0.780 kg
Cover box material:	Stainless steel
Application:	Incubating antibiotic test for raw milk



Temperature stability test results

1. Temperature stability test of LACTOSCAN TBS powered by 24V supply:

The Temperature stability test of LACTOSCAN TBS was made in Environmental chamber Astell, model EBC-006. The temperature of LACTOSCAN TBS was measured by using a calibrated digital thermometer -50 to 200 °C.

Environmental chamber's conditions:

- Temperature – 7 °C
- Humidity – 90%

LACTOSCAN TBS's conditions:


- Temperature - 7 °C
- Completely closed cover (using the no holes cover)
- Power supply 24V

LACTOSCAN TBS's operation period: 180 minutes

Results:

The results from the Temperature stability test of LACTOSCAN TBS shows that the device needs 5 (five) minutes to reach 47.5 °C. Then the temperature is stable during the rest of the 180 minutes period. The deviation is only ± 0.3 °C.

The controlled thermometer is placed in the thermometer's hole of LACTOSCAN TBS's plate.



2. Temperature stability test of LACTOSCAN TBS powered by 12V supply:

The Temperature stability test of LACTOSCAN TBS was made in Environmental chamber Astell, model EBC-006. The temperature of LACTOSCAN TBS was measured by using a calibrated digital thermometer -50 to 200 °C.

Environmental chamber's conditions:

- Temperature – 7 °C
- Humidity – 90%

LACTOSCAN TBS's conditions:


- Temperature - 7 °C
- Completely closed cover (using the no holes cover)
- Power supply 12V

LACTOSCAN TBS's operation period: 180 minutes

Results:

The results from the Temperature stability test of LACTOSCAN TBS shows that the device needs 12 (twelve) minutes to reach 47.5 °C. Then the temperature is stable during the rest of the 180 minutes period. The deviation is only ± 0.1 °C.

The controlled thermometer is placed in the thermometer's hole of LACTOSCAN TBS's plate.



WARRANTY CARD

LACTOSCAN TBS

Guarantee period is 1 (one) year after purchasing date. Improper handling, transport and storage will invalidate the guarantee. Guarantee is void if warranty labels are removed.

Serial №

Date of purchase:

Password:

Distributor:

Distributor:

Signature:

Stamp:

Milkotronic Ltd. reserves the right to select the methods used to analyze a product to determine the validity of the guarantee. No warranty is granted for products beyond their expiration date.

Contacts

For additional information and technical support, visit our website www.lactoscan.com or contact us:

Headquarters and Main service:

4 Narodni buditeli str.

8900 Nova Zagora

Bulgaria

Tel/fax: + 359 457 670 82

office@lactoscan.com